
MESI

Simplifying Diagnostics

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Přehled

Upozornění

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COMPARISON

Comparison between the handheld **Doppler** probe and **MESI** arterial assessment devices

MESI ABI measurement devices vs. the handheld Doppler probe method.

	Doppler Probe	MESI ABPI MD	MESI mTABLET ABI
Measurement duration	30 min¹	min	1 min
Pre-measurement resting	10-20 min⁰	min	0 min
Measuring process	One extremity at a time	Simultaneous 3-cuff measurement with	Simultaneous 3-cuff measurement with SmartArm™ detection
Additional education	YES	NO	NO
Calculations	Manual	Automatic^A	Automatic
Pulse waveform review	NO	YES in all extremities	YES in all extremities
Measurement report (PDF and print-out)	NO	Automatic via PC	Automatic from the mTABLET or any web-enabled device
Colour screen	NO	YES	YES large tablet with intuitive interface
Automated TBI extension	NO	NO	YES through the MESI mTABLET system
Gel application	YES	NO	NO
Clothes removal	YES	NO	NO
Colour-coded cuffs for easier placement	NO	YES	YES Pre-shaped cuffs in multiple sizes
Detection of severe PAD with a clear alert	NO	YES	YES
Direct storage into patient EHR	NO	NO	YES MESI mRECORDS platform

MESI mTABLET TBI vs. the handheld Doppler probe method.

	Doppler Probe	MESI mTABLET TBI
Measurement duration	30 min	1 min
Measuring process	One extremity at a time	Automatic simultaneous, 1-step measurement
Calculations	Manual	Automatic
Pulse waveform review	NO	YES oscillation graphs, pulse waveforms and PPG pulse waveforms
Measurement report (PDF and print-out)	NO	Automatic from the mTABLET or any web-enabled device
Direct storage into patient EHR	NO	YES MESI mRECORDS platform
Automated ABI extension	NO	YES through the MESI mTABLET system
Gel application	YES	NO
Clothes removal	YES	NO
Additional charge for separate software	YES	NO

MESI mTABLET TBI revolutionizes the way a toe-brachial index measurement is performed, shortening the procedure from 30 minutes to one minute. With a one-time cuff placement and a 1-step, simultaneous measurement using the MESI mTABLET's intuitive interface, required additional training is minimized. The Adaptive LED PPG probes and multiple-size toe cuffs, along with temperature and skin thickness detector ensure the result is accurate. Automatic calculation with in-depth analysis of oscillation graphs and pulse waveforms in the arms, and PPG pulse waveforms in the toes provide an excellent insight into the patient's state, even in very low perfusion.

REFERENCE

- 1 Sandison AT. Degenerative vascular disease in the Egyptian mummy. *Medical History*, 1962; 6:77-81.
- 2 Harvey W. De Motu Cordis et Sanguinis in Animalibus. *Oxford: Blackwell*, 1628.
- 3 Harvey W. Circulatione Sanguinis. *Oxford: Blackwell*, 1649.
- 4 Bollinger A, Eckert J, Rüttimann B, Becker F (2000). The “galloping” history of intermittent claudication. *Vasa* 29(4):295–299.
- 5 Brodie B. Lectures illustrative of various aspects of pathology and surgery. *London: Longmans*, 1846
- 6 Ejrup B. Tonoscillography after exercise in peripheral vascular disease and coarctation of the aorta. *Am Heart J* 1948; 35: 41-57.
- 7 Writing Committee Members (2017). 2016 AHA/ACC guideline on the management of patients with lower extremity peripheral artery disease: Executive Summary: A report of the American college of cardiology/ American Heart Association task force on clinical practice guidelines. *Circulation*, 135(12), e686-e725.
- 8 Criqui MH, Vargas V, Denenberg JO, Ho E, Allison M, Langer RD, Gamst A, Bundens WP, Fronck A. Ethnicity and peripheral arterial disease: the San Diego Population Study. *Circulation*, 2005; 112:2703–2707.
- 9 Hiramoto JS, Katz R, Weisman S, Conte M. Gender-specific risk factors for peripheral artery disease in a voluntary screening population. *J Am Heart Assoc.* 2014;3(2):e000651.
- 10 Joosten MM, Pai JK, Bertoia ML, et al. Associations between conventional cardiovascular risk factors and risk of peripheral artery disease in men. *JAMA.* 2012;308(16):1660–1667.
- 11 Conen D, Everett BM, Kurth T, et al. Smoking, smoking cessation, [corrected] and risk for symptomatic peripheral artery disease in women: a cohort study [published correction appears in *Ann Intern Med.* 2011 Sep 20;155(6):408]. *Ann Intern Med.* 2011;154(11):719–726.

- 12** Steven P. Marso, William R. Hiatt. Peripheral Arterial Disease in Patients With Diabetes. *Journal of the American College of Cardiology* Mar 2006, 47 (5) 921-929.
- 13** Olin JW, Sealove BA. Peripheral artery disease: current insight into the disease and its diagnosis and management. *Mayo Clin Proc.* 2010;85(7):678–692.
- 14** Denis L. Clement, Marc L. De Buyzere and Daniel A. Duprez. Hypertension in Peripheral Arterial Disease. *Current Pharmaceutical Design* (2004) 10: 3615.
- 15** Hicks CW, Yang C, Ndumele CE, et al. Associations of Obesity With Incident Hospitalization Related to Peripheral Artery Disease and Critical Limb Ischemia in the ARIC Study. *J Am Heart Assoc.* 2018;7(16):e008644.
- 16** Hirsch AT, Criqui MH, Treat-Jacobson D, et al. Peripheral Arterial Disease Detection, Awareness, and Treatment in Primary Care. *JAMA.* 2001;286(11):1317–1324.
- 17** McGrae McDermott M, Greenland P, Liu K, et al. Leg Symptoms in Peripheral Arterial Disease: Associated Clinical Characteristics and Functional Impairment. *JAMA.* 2001;286(13):1599–1606.
- 18** M. Condorelli, G. Brevetti. Intermittent claudication: an historical perspective. *European Heart Journal Supplements*, Volume 4, Issue suppl_B, 1 March 2002, Pages B2–B7.
- 19** Mohler III ER. Peripheral Arterial Disease: Identification and Implications. *Arch Intern Med.* 2003;163(19):2306–2314.
- 20** Paul W. Wennberg. Approach to the Patient With Peripheral Arterial Disease. *Circulation.* 2013;128:2241–2250.
- 21** Spentzouris G, Labropoulos N. The evaluation of lower-extremity ulcers. *Semin Intervent Radiol.* 2009;26(4):286–295.
- 22** Swaminathan A, Vemulapalli S, Patel MR, Jones WS. Lower extremity amputation in peripheral artery disease: improving patient outcomes. *Vasc Health Risk Manag.* 2014;10:417–424. Published 2014 Jul 16.
- 23** Polonsky TS, Taillon LA, Sheth H, Min JK, Archer SL, Ward RP. The association between erectile dysfunction and peripheral arterial disease as determined by screening ankle-brachial index testing. *Atherosclerosis.* 2009;7:440–444.

- 24** [www.who.int/news-room/fact-sheets/detail/cardiovascular-diseases-\(c-vds\)](http://www.who.int/news-room/fact-sheets/detail/cardiovascular-diseases-(c-vds))
- 25** www.who.int/news-room/fact-sheets/detail/cancer
- 26** Fowkes FG, Rudan D, Rudan I, Aboyans V, Denenberg JO, McDermott MM, Norman PE, Sampson UK, Williams LJ, Mensah GA, Criqui MH. Comparison of global estimates of prevalence and risk factors for peripheral artery disease in 2000 and 2010: a systematic review and analysis. *Lancet*. 2013 Oct 19;382(9901):1329-40.
- 27** www.who.int/healthinfo/global_burden_disease/metrics_daly/en/
- 28** Fowkes FG, Aboyans V, Fowkes FJ, McDermott MM, Sampson UK, Criqui MH. Peripheral artery disease: epidemiology and global perspectives. *Nat Rev Cardiol*. 2017 Mar;14(3):156-170.
- 29** J.P.Pell. Impact of intermittent claudication on quality of life. *European Journal of Vascular and Endovascular Surgery*. Volume 9, Issue 4, May 1995, Pages 469-472.
- 30** Meijer WT, Hoes AW, Rutgers D, Bots ML, Hofman A, Grobbee DE. Peripheral arterial disease in the elderly: The Rotterdam Study. *Arterioscler Thromb Vasc Biol*. 1998 Feb;18(2):185-92.
- 31** F G R FOWKES, E HOUSLEY, E H H CAWOOD, C C A MACINTYRE, C V RUCKLEY, R J PRESCOTT. Edinburgh Artery Study: Prevalence of Asymptomatic and Symptomatic Peripheral Arterial Disease in the General Population. *International Journal of Epidemiology*, Volume 20, Issue 2, June 1991, Pages 384–392.
- 32** Shamma NW. Epidemiology, classification, and modifiable risk factors of peripheral arterial disease. *Vasc Health Risk Manag*. 2007;3(2):229–234.
- 33** Poredoš, P., & Jug, B. (2007). The Prevalence of Peripheral Arterial Disease in High Risk Subjects and Coronary or Cerebrovascular Patients. *Angiology*, 58(3), 309–315.
- 34** Dieter, R. S., Tomasson, J., Gudjonsson, T., Brown, R. L., Vitcenda, M., Einerson, J., & McBride, P. E. (2003). Lower extremity peripheral arterial disease in hospitalized patients with coronary artery disease. *Vascular Medicine*, 8(4), 233–236.

- 35** Atmer B, Jogestrand T, Laska J, Lund F. Peripheral artery disease in patients with coronary artery disease. *Int Angiol.* 1995 Mar;14(1):89-93.
- 36** Cotter G, Cannon CP, McCabe CH, Michowitz Y, Kaluski E, Charlesworth A, Milo O, Bentley J, Blatt A, Krakover R, Zimlichman R, Reisin L, Marmor A, Lewis B, Vered Z, Caspi A, Braunwald E; OPUS-TIMI 16 Investigators. Prior peripheral arterial disease and cerebrovascular disease are independent predictors of adverse outcome in patients with acute coronary syndromes: are we doing enough? Results from the Orbofiban in Patients with Unstable Coronary Syndromes-Thrombolysis In Myocardial Infarction (OPUS-TIMI) 16 study. *Am Heart J.* 2003 Apr;145(4):622-7.
- 37** Cheng, S. W., Wu, L. L., Lau, H. , Ting, A. C. and Wong, J. (1999). PREVALENCE OF SIGNIFICANT CAROTID STENOSIS IN CHINESE PATIENTS WITH PERIPHERAL AND CORONARY ARTERY DISEASE. *Aust. N.Z. J. Surg.*, 69: 44-47.
- 38** Wilkins E, Wilson L, Wickramasinghe K, Bhatnagar P, Leal J, Luen-go-Fernandez R, Burns R, Rayner M, Townsend N (2017). European Cardiovascular Disease Statistics 2017. *European Heart Network*, Brussels.
- 39** Scully RE, Arnaoutakis DJ, DeBord Smith A, Semel M, Nguyen LL. Estimated annual health care expenditures in individuals with peripheral arterial disease. *J Vasc Surg.* 2018 Feb;67(2):558-567.
- 40** Pål Hasvold, Joakim Nordanstig, Björn Kragsterman, Thomas Kristensen, Mårten Falkenberg, Saga Johansson, Marcus Thuresson, Birgitta Sigvant. Long-term cardiovascular outcome, use of resources, and healthcare costs in patients with peripheral artery disease: results from a nationwide Swedish study. *European Heart Journal - Quality of Care and Clinical Outcomes*, Volume 4, Issue 1, January 2018, Pages 10–17.
- 41** Spentzouris G, Labropoulos N. The evaluation of lower-extremity ulcers. *Semin Intervent Radiol.* 2009;26(4):286–295.
- 42** Samuel R. Nussbaum, Marissa J. Carter, PhD, Caroline E. Fife, Joan DaVanzo, PhD, Randall Haight, Marcia Nusgart, Donna Cartwright. An Economic Evaluation of the Impact, Cost, and Medicare Policy Implications of Chronic Nonhealing Wounds. *Value in Health*, Volume 21, Issue 1, 27 – 32.

- 43** Dick F, Diehm N, Galimanis A, Husmann M, Schmidli J, Baumgartner I. Surgical or endovascular revascularization in patients with critical limb ischemia: influence of diabetes mellitus on clinical outcome. *J Vasc Surg*. 2007 Apr;45(4):751-61.
- 44** Adam DJ, Beard JD, Cleveland T, Bell J, Bradbury AW, Forbes JF, Fowkes FG, Gillespie I, Ruckley CV, Raab G, Storkey H; BASIL trial participants. Bypass versus angioplasty in severe ischaemia of the leg (BASIL): multicentre, randomised controlled trial. *Lancet*. 2005 Dec 3;366(9501):1925-34.
- 45** Stoyioglou A, Jaff MR. Medical treatment of peripheral arterial disease: a comprehensive review. *J Vasc Interv Radiol*. 2004 Nov;15(11):1197-207.
- 46** Abu Dabrh AM, Steffen MW, Undavalli C, Asi N, Wang Z, Elamin MB, Conte MS, Murad MH. The natural history of untreated severe or critical limb ischemia. *J Vasc Surg*. 2015 Dec;62(6):1642-51.e3.
- 47** Norgren L, Hiatt WR, Dormandy JA, Nehler MR, Harris KA, Fowkes FG; TASC II Working Group. Inter-Society Consensus for the Management of Peripheral Arterial Disease (TASC II). *J Vasc Surg*. 2007 Jan;45 Suppl S:S5-67.
- 48** Becker F1, Robert-Ebadi H, Ricco JB, Setacci C, Cao P, de Donato G, Eckstein HH, De Rango P, Diehm N, Schmidli J, Teraa M, Moll FL, Dick F, Davies AH, Lepäntalo M, Apelqvist J. Chapter I: Definitions, epidemiology, clinical presentation and prognosis. *Eur J Vasc Endovasc Surg*. 2011 Dec;42 Suppl 2:S4-12.
- 49** Barshes NR, Chambers JD, Cohen J, Belkin M. Model To Optimize Healthcare Value in Ischemic Extremities 1 (MOVIE) Study Collaborators. Cost-effectiveness in the contemporary management of critical limb ischemia with tissue loss. *J Vasc Surg*. 2012 Oct;56(4):1015-24.e1.
- 50** Alan T. Hirsch, Ziv J. Haskal, Norman R. Hertzler, Curtis W. Bakal, Mark A. Creager, Jonathan L. Halperin, Loren F. Hiratzka, William R.C. Murphy, Jeffrey W. Olin, Jules B. Puschett, Kenneth A. Rosenfield, David Sacks, James C. Stanley, Lloyd M. Taylor, Christopher J. White, John White, Rodney A. White, Elliott M. Antman, Sidney C. Smith, Cynthia D. Adams, Jeffrey L. Anderson, David P. Faxon, Valentin Fuster, Raymond J. Gibbons, Jonathan L. Halperin, Loren F. Hiratzka, Sharon A. Hunt, Alice K. Jacobs, Rick Nishimura, Joseph P. Ornato, Richard L. Page, Barbara Riegel. ACC/AHA 2005 Guidelines for the Management of Patients With Peripheral

Arterial Disease (Lower Extremity, Renal, Mesenteric, and Abdominal Aortic): A Collaborative Report from the American Association for Vascular Surgery/Society for Vascular Surgery, Society for Cardiovascular Angiography and Interventions, Society for Vascular Medicine and Biology, Society of Interventional Radiology, and the ACC/AHA Task Force on Practice Guidelines (Writing Committee to Develop Guidelines for the Management of Patients With Peripheral Arterial Disease). *Journal of the American College of Cardiology*. Mar 2006, 47 (6) e1-e192.

51 Lin, P. H., Bechara, C., Kougias, P., Huynh, T. T., LeMaire, S. A., & Coselli, J. S. (2009). Assessment of Aortic Pathology and Peripheral Arterial Disease Using Multidetector Computed Tomographic Angiography. *Vascular and Endovascular Surgery*. 42(6), 583–598.

52 O'Hare AM. Management of peripheral arterial disease in chronic kidney disease. *Cardiol Clin* 2005;23(3):225–236.

53 Brearley S, Shearman CP, Simms MH. *Peripheral pulse palpation: an unreliable physical sign*. *Ann R Coll Surg Engl*. 1992;74(3):169–171.

54 Kaiser, V & Kester, A.D.M. & Stoffers, Henri & J Kitslaar, P & Knottnerus, J.A.. (1999). The Influence of Experience on the Reproducibility of the Ankle–brachial Systolic Pressure Ratio in Peripheral Arterial Occlusive Disease. *European journal of vascular and endovascular surgery:the official journal of the European Society for Vascular Surgery*. 18. 25-9. 10.1053/ejvs.1999.0843.

55 Varetto G, Magnoni F, Aluigi L, Antignani PL, Ardita G, et al. (2019). Comparison of Ankle-Brachial Index (ABI) Measurement between a New Oscillometric Device (MESI ABPI Md®) and the Standard Doppler Method in the Diagnosis of Lower Extremity Arterial Disease (LEAD). *J Non Invasive Vasc Invest* 4: 012.

56 Aboyans et al., 2017 ESC Guidelines on the Diagnosis and Treatment of Peripheral Arterial Diseases, in collaboration with the European Society for Vascular Surgery (ESVS), chapter 10.2 Diagnostic tests.

57 María Martín-Fuentes, Lucrecia Herranz, Lourdes Saez-de-Ibarra, Luis Felipe Pallardo: Low Toe-Brachial Pressure Index Predicts Increased Risk of Recurrent Cardiovascular Disease in Type 2 Diabetes. *Diabetes Care* Apr 2015, 38 (4) e53-e54; DOI: 10.2337/dc14-1843.

58 Victor Aboyans, Jean-Baptiste Ricco, Marie-Louise E L Bartelink, Martin Björck, Marianne Brodmann, Tina Cohnert, Jean-Philippe Collet, Martin Czerny, Marco De Carlo, Sebastian Debus, Christine Espinola-Klein, Thomas Kahan, Serge Kownator, Lucia Mazzolai, A Ross Naylor, Marco Roffi, Joachim Röther, Muriel Sprynger, Michal Tendera, Gunnar Tepe, Maarit Venermo, Charalambos Vlachopoulos, Ileana Desormais, ESC Scientific Document Group. 2017 ESC Guidelines on the Diagnosis and Treatment of Peripheral Arterial Diseases, in collaboration with the European Society for Vascular Surgery (ESVS): Document covering atherosclerotic disease of extracranial carotid and vertebral, mesenteric, renal, upper and lower extremity arteries. The Task Force for the Diagnosis and Treatment of Peripheral Arterial Diseases of the European Society of Cardiology (ESC) and of the European Society for Vascular Surgery (ESVS). *European Heart Journal*, Volume 39, Issue 9, 01 March 2018, Pages 763–816.

Get your ABPI MD, full training
and local support from:

Promedeus

Promedeus

Maříkova 1899/1
621 00 Brno

objednavky@promedeus.cz
+420 702 248 624

MESI, Ltd.
Leskoškova cesta 11a
1000 Ljubljana, Slovenia, EU

www.mesimedical.com

E: info@mesimedical.com
T: +386 1 620 34 87

Facebook / **MESIdoo**
Twitter / **MESImedical**
Linkedin / **MESI**

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