



Ankalaps[®]

Easy Click System
Reusable Laparoscopic Instruments

T-Fix Retractor

Table Fixed Retractor Systems

our effort is for human...


www.ankasurgical.com



Ankalaps[®]



- Ergonomic
- Durable
- Precision
- Secure
- Easy to use



Easy Click System

Reusable Laparoscopic Instruments



T-Fix Retractor

Table Fixed Retractor Systems

Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Eseşehir Mah. Haseki Sok. Anka Plaza No:36 Ümraniye/İstanbul/Türkiye

Tel:+90 216 3150477 Fax:+90 216 3150478 Web:www.ankasurgical.com e-mail:info@ankasurgical.com



Table of Contents

Reusable Laparoscopic Instruments

Handles	Page 3
Dissectors 5-10mm	Page 4
Scissors 5mm	Page 4
Graspers 5mm	Page 4-5
Graspers 10mm	Page 6
Trocars	Page 7
Needle Holder	Page 8
Suction / Irrigation	Page 9
Monopolar HF Electrode	Page 9
Veress Needle	Page 10
Knot Pusher	Page 10
Fascia Closure Forceps	Page 10
Fan Retractors	Page 11
Laparoscopes	Page 11
Laparoscopic Holder and Tables	Page 25

Table Fixed Retractor Systems

T-Fix Retractor System / Nathanson

Model 1	Page 12
Model 2	Page 12

T-Fix Retractor System / Bookwalter

Clamps	Page 13-14
Horizontal Bar	Page 14
Rings	Page 14-15
Balfour Blades	Page 16
Malleable Blades	Page 16
Harrington Blades	Page 16
Kelly Blades	Page 17
Daever Blades	Page 17
Deep Pelvic Blades	Page 17
Adapters	Page 15-17

T-Fix Retractor System / Thompson

Table Connection Arm	Page 18
Bar Clamp	Page 19
Crossbars	Page 19-20
Blade Holders	Page 20
Balfour Blades	Page 21
Harrington Blades	Page 21
Fenestrated Blades	Page 21
Kelly Blades	Page 21
Splanchnic Blades	Page 22
Fence Blades	Page 22
Malleable Blades	Page 22
Malleable Finger Blades	Page 23
Malleable Slotted Blades	Page 23
Double Angled Blades	Page 23
Heavy Blades	Page 23
ST Marks Blades	Page 24
Deaver Blades	Page 24
Meyerding Blades	Page 24
Notes	Page 26



Ankalaps[®]

Handpiece 5mm

- Monopolar cautery pin
- Easy ratchet for right or left handed operation
- Easy Disassembly Button
- Comfortable single-handed rotation knob
- Washing Port for fast, safe, and effective cleaning
- High dielectric strength shaft material



Precision Performance

- The reusable system offers a superior experience for surgeons and offers an environmentally friendly alternative to single-use instruments. It is comfortable and provides tactile feedback.
- Available in shaft lengths of 21, 33 and 42cm to accommodate patient anatomy and surgical procedures.
- Maximum precision, intuitive control for clinical, ergonomic design that can be divided into 3 parts

Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Esenşehir Mah. Haseki Sok. Anka Plaza No:36 Ümraniye/İstanbul/Türkiye

Tel:+90 216 3150477 Fax:+90 216 3150478 Web:www.ankasurgical.com e-mail:info@ankasurgical.com



Dissectors



101.011B
Right Angle Dissector 5mm
Working Tip Length is 17.0mm
True right angle for wide tissue spread and dissection around delicate structures



101.013B
Maryland Dissector 5mm
Working Tip Length is 20.0mm
A slender jaw design with micro-serrations and a pronounced angle at the tip allows for more delicate tissue manipulation



101.011A
Right Angle Dissector 10mm
Working Tip Length is 17.0mm
Extra tip length for greater visibility. Requires a 10mm Trocar.

Scissors



101.014
Metzenbaum Scissor 10mm
Curved Tip Length is 17.0mm



101.017
Straight Scissor 5mm
Length is 17.0mm



101.019
Metzenbaum Scissor 5mm
Curved Tip Length is 17.0mm

Graspers



101.035
Babcock Grasper 5mm
Working Tip Length is 35.0mm
Traditional pattern with larger grasping pads to permit more secure, gentle grasping of delicate tissues.



101.036
Fenestrated Grasper 5mm
Working Tip Length is 22.0mm
Fenestrated and grooved slot to allow for strong yet gentle grasping.



101.040B
Long Bowel Grasper 5mm
Working Tip Length is 36.0mm
Longer jaws promote atraumatic grasping and manipulation of bowel, stomach and mesentery. Micro-serrations.



101.040B2
Long Fenestrated Grasper 5mm
Working Tip Length is 36.0mm
Longer jaws promote atraumatic grasping and manipulation of bowel, stomach and mesentery. Slotted with micro-serrations.



101.041A
Lapclinch Grasper 5mm
Working Length 22.0 mm
End design of large gear. Strong grip in thick tissues



101.041C
Modified Lapclinch Grasper 5mm
Working Length 18.0 mm
End design of large gear. Strong grip in thick tissues



Graspers



101.043A

Cobra Tooth Grasper 5mm

Working Type Length is 17.0mm
Jaw design provides an extremely secure grip. Improved traditional design



101.043D

Cobra Tooth Grasper 5mm

Working Type Length is 12.0mm
Jaw design provides an extremely secure grip. Improved traditional design Single action.



101.046

Max. Atravmatic Grasper 5mm

Working Type Length is 17.0mm
Precision and robust grip



101.047A

Reversing Tooth Grasper 5mm

Working Type Length is 22.0mm



101.054

Curved Dissecting Knot Grasper 5mm

Working Type Length is 17.0mm



101.100

Claw Grasper 5mm

Working Type Length is 22.0mm



101.102A

Stomach Grasper 5mm

Working Type Length is 22.0mm



101.102B

Stomach Grasper 5mm

Working Type Length is 22.0mm



101.102C

Stomach Grasper 5mm

Working Type Length is 22.0mm



101.108

Fixed Grasper 5mm

Working Type Length is 22.0mm

Reliability



Graspers 10mm



101.006B
Babcock Grasper 10mm
Double Action



101.006
Claw Grasper 10mm
Single Action



101.007A
Atravmatic Grasper 10mm
Double Action



101.079
Anvil Grasper 10mm
Single action



201.004B
Single Tooth Tenaculum Forceps
5mm X 420mm

201.003B
Single Tooth Tenaculum Forceps
10mm X 420mm

Precision Performance

Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Esenşehir Mah. Haseki Sok. Anka Plaza No:36 Ümraniye/İstanbul/Türkiye

Tel:+90 216 3150477 Fax:+90 216 3150478 Web:www.ankasurgical.com e-mail:info@ankasurgical.com



Trocars


101.008A

Trocar Point Head- Smooth Cannula
10,5mmx100mm - 10,5mmx120mm



101.008B

Trocar Point Head- Smooth Cannula
5,5mmx100mm - 5,5mmx120mm



Magnetic Sheet Valve

101.008C

Protected head- Smooth Cannula
10,5mmx100mm - 10,5mmx120mm



101.008D

Protected head- Smooth Cannula
5,5mmx100mm - 5,5mmx120mm



101.071A

Reducer
Ø10mm/Ø 5 mm



101.008E

Protected head- Screw Cannula
10,5mmx100mm - 10,5mmx120mm



101.008F

Protected head- Screw Cannula
5,5mmx100mm - 5,5mmx120mm



101.071D

Reducer
Ø10mm/Ø 5 mm





Needle Holder



101.023A

Straight

5x330mm - 5x420mm



101.023B

Curved

5x330mm - 5x420mm



101.023C

Self Righting

5x330mm - 5x420mm



101.023D

Self Righting

5x330mm - 5x420mm



101.023F

Straight

5x330mm - 5x420mm



101.023E

Curved

5x330mm - 5x420mm



101.023G

Self Righting

5x330mm - 5x420mm



101.023H

Self Righting

5x330mm - 5x420mm



101.023I

Self Righting

5x330mm - 5x420mm



101.023J

Straight

5x330mm - 5x420mm



101.023L

Curved

5x330mm - 5x420mm



101.023K

Self Righting

5x330mm - 5x420mm



101.023M

Self Righting

5x330mm - 5x420mm



101.023N

Self Righting

5x330mm - 5x420mm



101.025A

Straight

5x330mm - 5x420mm



101.025B

Curved

5x330mm - 5x420mm



101.025C

Self Righting

5x330mm - 5x420mm



101.025D

Righting

5x330mm - 5x420mm



101.026A

Straight

5x330mm - 5x420mm



101.026B

Curved

5x330mm - 5x420mm



101.025C

Self Righting

5x330mm - 5x420mm



101.025D

Self Righting

5x330mm - 5x420mm



101.025E

Self Righting

5x330mm - 5x420mm



Suction / Irrigation

101.068

Sliding Model Ø5mm- Ø10mm
340mm - 420mm



101.069

Trompet Valve Model Ø5mm- Ø10mm
340mm - 420mm



Suction Tube Ø5mm
340mm - 420mm



Suction Tube with aspiration needle Ø5mm
340mm - 420mm



Suction Tube Ø10 mm
340mm - 420mm



101.083B

Sliding Model with Hook Probe
Ø5mm 340mm - 420mm



Monopolar HF Electrode



101.027
Hook Electrode



101.027D
Hook Electrode



101.027F
Hook Electrode



101.027E
Hook Electrode



101.030
Spatula Electrode



101.031
Needle Electrode



101.028
Ballpoint Electrode



101.029
Knife Electrode



101.062B
Monopolar Cable


Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Eseneşhir Mah. Haseki Sok. Anka Plaza No:36 Ümraniye/İstanbul/Türkiye

Tel:+90 216 3150477 Fax:+90 216 3150478 Web:www.ankasurgical.com e-mail:info@ankasurgical.com



Veress Needle




Fascia Closure Forceps



102.060B
Up shed
 $\varnothing 2,4 \times 120$ mm

102.060C
Under shed
 $\varnothing 2,4 \times 120$ mm






Knot Pusher






Fan Retractors



Laparoscopes

AUTOCLAVABLE

KE-01A 30° / Ø10mm x 330mm



KE-01B 0° / Ø5mm x 315mm






Nathanson

Surgical Holders, Positioning Arms, and Laparoscopic Retractors

Model 1



REF	PART	SIZES	WIDTH	ANGLE
NE-101	HOOK BLADE	105 mm	110 mm	75°
NE-102	HOOK BLADE	165 mm	110 mm	90°
NE-103	HOOK BLADE	210 mm	185 mm	110°
NE-104	HOOK BLADE	220 mm	110 mm	125°
NE-105	HOOK BLADE	235 mm	185 mm	135°
NE-106	HOOK HOLDER	500 mm		
NE-107	TABLE FIXING ROD	650 mm		
NE-108	TABLE FIXING ROD ADAPTOR	40x12x20 mm		
NE-109	OPTICAL HOLDER	1000 mm		
NE-110	STERILISATION CONTAINER	580x280x260 mm		
NE-100	COMPLETE SET			

**Faster Setup,
Fewer Parts,
Less Weight,
More Adjustable**



Model 2




REF	PART	SIZES	ANGLE
NE-101 P	HOOK BLADE 1	105 mm	75°
NE-102 P	HOOK BLADE 2	165 mm	90°
NE-103 P	HOOK BLADE 3	210 mm	110°
NE-104 P	HOOK BLADE 4	220 mm	125°
NE-105 P	HOOK BLADE 5	235 mm	135°
NE-106 P	HOOK HOLDER	500 mm	
NE-107 P	TABLE FIXING ROD	650 mm	
NE-108 P	TABLE FIXING ROD ADAPTOR	40x12x20mm	
NE-109 P	OPTICAL HOLDER	500 mm	
NE-110 P	STERILISATION CONTAINER	580x280x200 mm	
NE-100 P	COMPLETE SET		





Bookwalter

**-Robust,
-Modular,
-Ergonomic**



Ankalaps T-Fix Retractor Bookwalter model is manufactured from 316 quality stainless steel by our company. Our company guarantees for 5 years against all kinds of manufacturing defects. Parts manufactured in standard are listed below. Parts can be created by selecting the required set of retractors. In addition to standard set parts, optional sizes can be produced. As with all surgical instruments, it is absolutely necessary to make annual maintenance of the products in certain periods. Each generated set has a unique code number. For each set, maintenance activities are carried out free of charge once a year for 5 years.

Flat Table Clamp

Item No	Overall Length
BW-101-460	460 mm
BW-101-610	610 mm



Blade Clamp

Item No	Description
BW-102-001	Non Tilting
BW-102-002	Tilting





Bookwalter


Post Clamp

Item No	Description
BW-102-003	For 25mm Bar



Horizontal Bar

Item No	Length
BW-103A-480	480 mm
BW-103A-560	560 mm
BW-103B-380	380 mm
BW-103B-560	560 mm



Oval Ring

Item No	Description	Width/Length
BW-104A-203-355	Small	215 x 406mm
BW-104A-265-380	Medium	265 x 380 mm
BW-104A-315-470	Large	315 x 470 mm



Segmented Oval Ring

Item No	Description	Width/Length
BW-104B-203-355	Small	203 x 355mm
BW-104B-255-495	Medium	255 x 495 mm
BW-104B-305-495	Large	305 x 495 mm
BW-104B-406-597	X-Large	406 x 597 mm





Bookwalter

Round Ring

Item No	Description	Diameter
BW-105A-215	Small	215 mm
BW-105A-267	Medium	267 mm
BW-105A-317	Large	317 mm



Segmented Round Ring

Item No	Description	Diameter
BW-105B-215	Small	215 mm
BW-105B-267	Medium	267 mm
BW-105B-317	Large	317 mm
BW-105B-406	X-Large	406 mm



Segmented Oval and Round Ring Pair

Item No	Lenght
BW-105-102	102 mm
BW-105-152	152 mm
BW-105-203	203 mm



Segmented Ring Screw

Item No	Description
BW-105-001	Replacement screw for T-Fix Segmented Rings



-Robust,
-Modular,
-Ergonomic





Bookwalter

Balfour Blades

Item No	Width / Length
BW-106-76-63	76mm x 63mm
BW-106-76-83	76mm x 83mm
BW-106-76-90	76mm x 90mm



Malleable Blades

Item No	Item No	Width / Length
Angled	Flat	
BW-107-20-51	BW-107B-20-51	20mm x 51mm
BW-107-20-76	BW-107B-20-76	20mm x 76mm
BW-107A-26-26	BW-107B-26-26	26mm x 26mm
BW-107A-26-102	BW-107B-26-102	26mm x 102mm
BW-107A-26-127	BW-107B-20-127	20mm x 127mm
BW-107A-38-76	BW-107B-38-76	38mm x 76mm
BW-107A-38-102	BW-107B-38-102	38mm x 102mm
BW-107A-38-152	BW-107B-38-152	38mm x 152mm
BW-107A-38-203	BW-107B-38-203	38mm x 203mm
BW-107A-51-152	BW-107B-51-152	51mm x 152mm
BW-107A-51-203	BW-107B-51-203	51mm x 203mm
BW-107A-51-254	BW-107B-51-254	51mm x 254mm
BW-107A-76-152	BW-107B-76-152	76mm x 152mm
BW-107A-76-203	BW-107B-76-203	76mm x 203mm
BW-107A-76-254	BW-107B-76-254	76mm x 254mm
BW-107A-76-305	BW-107B-76-305	76mm x 305mm



Harrington Blades

Item No	Width / Length
BW-108-38-127	38mm x 127mm
BW-108-63-140	63mm x 140mm
BW-108-89-254	89mm x 254mm
BW-108-114-254	114mm x 254mm





Bookwalter

Kelly Blades

Item No	Width / Length
BW-109-25-38	25mm x 38mm
BW-109-25-51	25mm x 51mm
BW-109-51-38	51mm x 38mm
BW-109-51-63	51mm x 63mm
BW-109-51-76	51mm x 76mm
BW-109-51-102	51mm x 102mm
BW-109-51-127	51mm x 127mm
BW-109-51-152	51mm x 152mm
BW-109-51-203	51mm x 203mm



Daever Blades

Item No	Width / Length
BW-110-25-114	25mm x 114mm
BW-110-38-114	38mm x 114mm
BW-110-51-114	51mm x 114mm



Deep Pelvic Blade

Item No	Width / Length
BW-111-60-190	60mm x 190mm




Adapter for Hand Held Retractors

Item No	Width / Length
BW-112	





Thompson



- General Surgery
- Bariatric Surgery
- Colon / Rectal Surgery
- Liver Surgery
- Oncologic Surgery
- Vascular Surgery
- Pediatric Surgery
- Retractor Set Systems for Small Incision Surgery

Ankalaps T-Fix Retractor Thomson model is manufactured from 316 quality stainless steel by our company. Our company guarantees for 5 years against all kinds of manufacturing defects. Parts manufactured in standard are listed below. Parts can be created by selecting the required set of retractors. In addition to standard set parts, optional sizes can be produced. As with all surgical instruments, it is absolutely necessary to make annual maintenance of the products in certain periods. Each generated set has a unique code number. For each set, maintenance activities are carried out free of charge once a year for 5 years.

Table Connection Arm

Item No	Description
TH-101-20-250	25 cm
TH-101-20-410	41 cm
TH-101-20-450	45 cm
TH-101-20-550	55 cm
TH-101-20-650	65 cm






Thompson

Bar Clamp

Item No	Description
TH-102-001	Bar Clamp 1 12mm x 12mm
TH-102-002	Bar Clamp 2 12mm x 12mm




Bilateral Crossbar

Item No	Description
TH-103-001	Bilateral Crossbar 184 x 343 x 184 mm
TH-103-002	Bilateral Crossbar Hinged 292 x 216 x 292 mm
TH-103-003	Bilateral Crossbar Hinged 368 x 254 x 368 mm



Unilateral Crossbar

Item No	Description
TH-104-150-500	Unilateral Crossbar 1 15x50cm
TH-104-280-330	Unilateral Crossbar 2 28x33cm
TH-104-200-400	Unilateral Crossbar 3 20x40cm
TH-104-250-250	Unilateral Crossbar 4 25x25cm
TH-104-200-250	Unilateral Crossbar 5 20x25cm
TH-104-250-180	Unilateral Crossbar 6 25x18cm
TH-104-150-150	Unilateral Crossbar 7 15x15cm






Thompson


Straight Crossbar

Item No	Description
TH-105-300	Straight Crossbar 1 30cm
TH-105-500	Straight Crossbar 2 50cm




Blade Holder Fixed

Item No	Description
TH-106-001F	Blade Holder Fixed



Blade Holder with Adjustable Angle

Item No	Description
TH-106-002A	ADJUSTABLE ANGLE






Thompson Blades

Balfour Blades

Item No	Item No	Description
Non-Lips	Lips	Non-Lips & Lips
TH-107A-32-32	TH-107B-32-32	32 x 32 mm
TH-107A-51-54	TH-107B-51-54	51 x 54 mm
TH-107A-76-51	TH-107B-76-51	76 x 51 mm
TH-107A-70-73	TH-107B-70-73	70 x 73 mm
TH-107A-83-73	TH-107B-83-73	83 x 73 mm
TH-107A-83-123	TH-107B-83-123	83 x 123 mm
TH-107A-83-175	TH-107B-83-175	83 x 175 mm



Harrington (Sweetheart) Blades

Item No	Description
TH-108H-38-127	38 x 127 mm
TH-108H-64-102	64 x 102 mm
TH-108H-64-152	64 x 152 mm
TH-108H-64-203	64 x 203 mm
TH-108H-64-254	64 x 254 mm
Malleable	Malleable
TH-108M-38-203	38 x 203 mm
TH-108M-64-152	64 x 152 mm
TH-108M-114-203	114 x 203 mm



Fenestrated Blades

Item No	Description
TH-109-44-89	44 x 89 mm
TH-109-48-27	48 x 27 mm



Kelly Blades

Item No	Description
TH-110-38-51	38 x 51 mm
TH-110-51-64	51 x 64 mm
TH-110-64-76	64 x 76 mm
TH-110-76-89	76 x 89 mm



Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Esenşehir Mah. Haseki Sok. Anka Plaza No:36 Ümraniye/İstanbul/Türkiye

Tel:+90 216 3150477 Fax:+90 216 3150478 Web:www.ankasurgical.com e-mail:info@ankasurgical.com



Thompson Blades

Splanchnich Blades

Item No	Description
TH-111-25-152	25 x 152 mm
TH-111-25-203	25 x 203 mm
TH-111-51-114	51 x 114 mm
TH-111-51-152	51 x 152 mm
TH-111-51-203	51 x 203 mm



Fence Blades

Item No	Description
TH-112-102-102	102 x 102 mm
TH-112-102-127	102 x 127 mm
TH-112-102-152	102 x 152 mm
TH-112-102-203	102 x 203 mm



Malleable Blades

Item No	Description
TH-113-13-90	13 x 90 mm
TH-113-13-64	13 x 64 mm
TH-113-13-152	13 x 152 mm
TH-113-13-203	13 x 203 mm
TH-113-25-104	25 x 104 mm
TH-113-25-110	25 x 110 mm
TH-113-25-127	25 x 127 mm
TH-113-25-254	25 x 254 mm
TH-113-38-203	38 x 203 mm
TH-113-38-254	38 x 254 mm
TH-113-51-102	51 x 102 mm
TH-113-51-203	51 x 203 mm
TH-113-51-254	51 x 254 mm
TH-113-64-254	64 x 254 mm
TH-113-76-254	76 x 254 mm
TH-113-102-25	102 x 254 mm






Thompson Blades

Malleable Finger Blades

Item No	Description
TH-114-76	76 mm
TH-114-102	102 mm
TH-114-152	152 mm
TH-114-203	203 mm



Malleable Slotted Blades

Item No	Description
TH-115-51-203	51 x 203 mm
TH-115-76-203	76 x 203 mm
TH-115-51-165	51 x 165 mm
TH-115-76-165	76 x 165 mm



Double Angled Blades

Item No	Description
TH-116-13-64	13 x 64 mm
TH-116-13-70	13 x 70 mm
TH-116-13-76	13 x 76 mm
TH-116-25-64	25 x 64 mm
TH-116-25-76	25 x 76 mm



Heany Blades

Item No	Description
TH-117-25-102	25 x 102 mm
TH-117-25-127	25 x 127mm





Thompson Blades

ST Marks Blades

Item No	Description
TH-118-64-140	64 x 140 mm (no lip)
TH-118-64-178	64 x 178 mm (no lip)
TH-118L-64-140	64 x 140 mm (with lip)
TH-118L-64-178	64 x 178 mm (with lip)



Deaver Blades (Non Malleable)

Item No	Description
TH-119-25-127	25 x 127 mm
TH-119-38-127	38 x 127 mm
TH-119-51-127	51 x 127 mm
TH-119-64-127	64 x 127 mm
TH-119-76-127	76 x 127 mm
TH-119-102-127	102 x 127 mm



Meyerding Blades

Item No	Description
TH-120-16-55	16 x 55 mm
TH-120-16-75	16 x 75 mm
TH-120-16-100	16 x 100 mm
TH-120-25-60	25 x 60 mm
TH-120-25-75	25 x 75 mm
TH-120-25-100	25 x 100 mm
TH-120-25-125	25 x 125 mm
TH-120-51-60	51 x 60 mm
TH-120-51-75	51 x 75 mm
TH-120-51-100	51 x 100 mm
TH-120-51-125	51 x 125 mm



Warranty Information

Ankalaps instruments are covered by one year warranty for defects in materials and workmanship. Items returned by Anka Sağlık Gereçleri for warranty service which are found to be defective in material or workmanship shall be repaired or replaced free of charge. This warranty is invalid if used for non-intended purposes, misused or accidentally damaged or used by unauthorized persons.





EasyClick Laparoscopic Instrument Retention System




For use in laparoscopic procedures,

- Holds all standard laparoscopic hand tools,
- Silicone tool holders protect and stabilize tools,
- Provides quick and easy access to instruments wherever and whenever you need them,
- Specially designed to hold 13 or more laparoscopic handpieces of any size, including clip applicators, laparoscopic hand tools, needle holders and stapling devices.
- The console type rack is attached directly to the surgical table rail. It provides a fast and safe presentation of the devices and moves with the surgical table.
- Keeps the instruments sterile. Prevents disruption of procedure to replace tools. Provides immediate access to instruments without losing sight of the monitor. The surgeon continues to focus on the procedure.
- The surgeon chooses the instrument as he wishes.

EasyClick Laparoscopic Instrument Organizer System

Laparoscopic hand tools container systems,

- Ensures regular placement of all standard laparoscopic hand tools in sterilization containers.
- Hand tools within the container multiplication, nesting, etc. reasons to prevent damage.
- Silicone tool holders protect and stabilize the tools.
- Provides quick and easy access to instruments wherever and whenever you need them.
- Prevents the tools from falling or damaging during insertion or removal in the sterilization container.





Anka Sağlık Gereçleri San. İç ve Dış Tic. Ltd. Şti

Esenşehir Mah.Haseki Sok. Anka Plaza No:36 Ümraniye - İstanbul/TURKEY

Tel:+90 216 315 04 77 Fax: +90 216 315 04 78

web:www.ankasurgical.com / e-mail:info@ankasurgical.com

our effort is for human...